

Protocol for passaging of HBEC3-KT

Version: May 2021

Evercyte Ord. No.:	CkHT-004-0230
Designation:	HBEC3-KT, human bronchial epithelial progenitor cells
Growth medium:	Keratinocyte-SFM Kit (Gibco, Cat# 17005-042):

Final components:

Keratinocyte-SFM basal medium (Gibco, Cat# 17005042)

Bovine Pituitary Extract (Gibco, Cat# 17005042)

EGF (Gibco, Cat# 17005042)

- take one bottle of Keratinocyte-SFM (1x) basal medium (500 ml)
- add 25 mg of BPE (Keratinocyte-SFM component)
- add 2.5 µg of EGF (Keratinocyte-SFM component)

- mix properly,
- store at 4°C for a maximum of 4 weeks
- temper the medium to room temperature (not 37°C) before use

Coating:	0.1 % Gelatin solution from porcine skin – Type A
----------	---

The coating solution is prepared by mixing the following components:

Powdered Gelatin from porcine skin – Type A (Sigma-Aldrich, Cat# G1890, stored at 4°C)

Cell Culture Grade Water (HyClone, Cat# SH30529.03, ready-to-use, stored at RT)

- weigh 2 g of Gelatin in glass bottle
- add 200 ml Cell Culture Grade Water
- transfer bottle to water bath to dissolve Gelatin
- autoclave resulting 1 % Gelatin solution
- aliquot 5 ml each
- store at 4°C

For coating of a T25 roux flask proceed as follows:

- liquefy the Gelatin solution at 37°C
- add 45 ml of Cell Culture Grade Water to 5 ml 1 % Gelatin solution and mix carefully
- store resulting 0.1 % Gelatin solution at 37°C
- transfer 2 ml of diluted Gelatin solution (0.1 %) to the T25 roux flask (80 µl/cm²)
- completely wet the surface of the culture flask
- incubate the culture flask at 37°C at least for 4 hours (up to one week)

- remove excess of Gelatin solution
- use culture flask immediately for seeding of cells, the surface must not dry out

Additional reagents: 0.025 % Trypsin-EDTA (Gibco, Cat# R-001-100, ready-to-use, stored at 4°C)
Defined Trypsin Inhibitor (Gibco, Cat# R007100, ready-to-use, stored at 4°C)
PBS (Sigma-Aldrich, Cat# D8537, ready-to-use, stored at RT)
0.1 % Gelatin from porcine skin – Type A (Sigma-Aldrich, Cat# G1890)
2 % FCS in PBS (FBS: Sigma-Aldrich, Cat# F7524, stored at 4°C)

- Passaging of cells:
- the new culture flasks have to be coated with Gelatin as described above
 - remove and discard the culture medium
 - wash the cells twice with PBS (each 160 $\mu\text{l}/\text{cm}^2$), remove PBS completely
 - add Trypsin-EDTA solution (20 $\mu\text{l}/\text{cm}^2$), make sure that all cells have been in contact with this solution
 - incubate the culture flask at 37°C for approximately 2-3 min
 - observe cell detachment under an inverted microscope
 - as soon as all cells are detached (if necessary, agitate the cells by gently hitting the flask) halt the Trypsin action by addition of 2 % FCS in PBS (220 $\mu\text{l}/\text{cm}^2$)
 - centrifuge at 170 g for 5 min
 - discard the supernatant, resuspend the cell pellet in the remaining droplet and add growth medium (about 160 $\mu\text{l}/\text{cm}^2$)
 - transfer appropriate aliquots of the cell suspension to Gelatin-coated culture vessels supplemented with growth medium (final volume of 240 $\mu\text{l}/\text{cm}^2$)
 - a split ratio of 1:4 twice a week is recommended (after having reached about 80-90 % confluence)
 - perform a medium change after 2-3 days if cells have not reached required cell density
 - cultivate cells at 37°C in a humidified atmosphere with 5% CO₂
-