

## Protocol for passaging of HUVEC/TERT66

Version: May 2021

Evercyte Ord. No.:	CHT-006-0066
Designation:	HUVEC/TERT66, umbilical vein endothelial cells
Growth medium:	<p>Endopan 300 SL Kit (PAN Biotech) supplemented with Panexin SL-S (PAN Biotech) and G418:</p> <p><u>Final components:</u></p> <p>Endopan 300 SL (PAN Biotech, Cat# P04-0065K without GA)</p> <p>Serum substitute Panexin SL-S (PAN Biotech, Cat# P04-0065S)</p> <p>20 µg/ml G418 (InvivoGen, Cat# ant-gn-5, 100 mg/ml stock solution, ready-to-use)</p> <ul style="list-style-type: none"> <li>- take one bottle (500 ml) of Endopan 300 SL basal medium</li> <li>- add all supplements from Endopan 300 SL Kit but GA</li> <li>- add one additional vial of Panexin SL-S (25 mL)</li> <li>- add 100 µl of G418 stock solution</li> </ul> <ul style="list-style-type: none"> <li>- mix properly</li> <li>- store at 4°C for a maximum of 4 weeks</li> <li>- temper the medium to room temperature (not 37°C) before use</li> </ul>
Additional reagents:	<p>PBS (Sigma-Aldrich, Cat# D8537, ready-to-use, stored at RT)</p> <p>0.1 % Gelatin (Sigma-Aldrich, Cat# G1393, 2 %), diluted in PBS</p> <p>0.05 % Trypsin-EDTA (Gibco, Cat#25300-054, ready-to-use, stored at 4°C after thawing)</p> <p>Defined Trypsin Inhibitor (Gibco, Cat# R007100, ready-to-use, stored at 4°C after thawing)</p>
Coating:	<p>0.1 % Gelatin solution</p> <p>The coating solution is prepared by mixing the following components:</p> <p>Gelatin (Sigma-Aldrich, Cat# G1393, 2 %, stored at 4°C)</p> <p>PBS (Sigma-Aldrich, Cat# D8537, ready-to-use, stored at RT)</p> <ul style="list-style-type: none"> <li>- liquify the Gelatin stock solution (2 %) at 37°C</li> <li>- add 38 ml PBS to 2 ml Gelatin stock solution (2 %)</li> <li>- mix carefully, store diluted Gelatin solution (0.1 %) in aliquots at 37°C</li> </ul> <p>For coating of a T25 roux flask proceed as follows:</p> <ul style="list-style-type: none"> <li>- transfer 2 ml of Gelatin solution (0.1 %) to a T25 roux flask (final 80 µl/cm<sup>2</sup>)</li> <li>- completely wet the surface of the culture flask</li> <li>- incubate the culture flask at 37°C between 10-60 min</li> <li>- remove excess of Gelatin solution</li> </ul>

---

use culture flask immediately for seeding of cells, the surface must not dry out

---

- Passaging of cells:
- the new culture flasks have to be coated as described above
  - remove and discard the culture medium
  - wash the cells twice with PBS (each 160  $\mu\text{l}/\text{cm}^2$ ), remove PBS completely
  - add 0.05 % Trypsin-EDTA solution (20  $\mu\text{l}/\text{cm}^2$ ), make sure that all cells have been in contact with this solution
  - incubate the culture flask at 37°C for approximately 3 min
  - observe cell detachment under an inverted microscope
  - as soon as all cells are detached (if necessary, agitate the cells by gently hitting the flask) add Defined Trypsin Inhibitor (about 20  $\mu\text{l}/\text{cm}^2$ )
  - resuspend the cells in growth medium (about 160  $\mu\text{l}/\text{cm}^2$ ) and aspirate the cells by pipetting
  - centrifuge at 170 g for 5 min
  - discard the supernatant, resuspend the cell pellet in the remaining droplet and add growth medium
  - add appropriate aliquots of the cell suspension to Gelatin coated culture vessels supplemented with growth medium (final volume of 240  $\mu\text{l}/\text{cm}^2$ )
  - a split ratio of 1:2 twice a week is recommended (after having reached about 90-95 % confluence)
  - cultivate cells at 37°C in a humidified atmosphere with 5% CO<sub>2</sub>

---

Related products: HUVEC/TERT2, umbilical vein endothelial cells (Evercyte, Cat# CHT-006-0008)

---