

Protocol for preparation of ColoUp medium

page 1 of 1

Version: May 2021

| Evercyte Ord. No.: | MHT-039 |
|--------------------|---|
| Reagents: | The ColoUp medium for cultivation of HCEC-1CT cells can either be ordered from Evercyte as ready-to-use medium (Cat# MHT-039) or as basal medium (Cat# MHT-039-B) plus supplements (Cat# MHT-039-S). The medium can also be prepared by mixing the following components: |
| | DMEM (Gibco, Cat # 10566016) / M199 (Gibco, Cat# 31150022) (4+1) 2 % Cosmic Calf Serum (Hyclone, Cat# SH30087, ready-to-use) 20 ng/ml hEGF (Sigma-Aldrich, Cat# E9644) 10 μg/ml Insulin (Sigma-Aldrich, Cat# 19278, ready-to-use) 2 μg/ml Apo-Transferrin (Sigma-Aldrich, Cat# T2036) 5 nM Sodium-Selenite (Sigma-Aldrich, Cat# S5261) 1 μg/ml Hydrocortisone (Sigma-Aldrich, Cat# H0396) - take one bottle (500 ml) of DMEM and discard 100 ml - add 100 ml of M199 and mix properly - discard 11.6 ml DMEM/M199 mixture - add 10 ml Cosmic Calf Serum (ready-to-use) - add 500 μl of hEGF stock (20 μg/ml, prepared in cell culture grade water) - add 500 μl of Insulin (10 mg/ml, ready-to-use) - add 25 μl of Sodium-Selenite stock (10 mg/ml, prepared in cell culture grade water) - add 500 μl Hydrocortisone stock (1 mg/ml, prepared in cell culture grade water) - mix properly and store at 4°C for up to 1 month - temper the medium to room temperature (not 37°C) before use |
| Related products: | ColoUp ready-to-use medium, 500 ml (Evercyte, Cat#MHT-039) ColoUp basal medium, 500 ml (Evercyte, Cat# MHT-039-B) ColoUp supplements (Evercyte, Cat# MHT-039-S) HCEC-1CT, human colon epithelial cells (Evercyte, Cat# CkHT-039-0229) |