

Protocol for cryopreservation of hTEC/SVTERT24-B

Version: May 2021

| | |
|----------------------|--|
| Evercyte Ord. No.: | CLHT-010-0024-B |
| Designation: | hTEC/SVTERT24-B, human thymic epithelial cells |
| Freezing medium: | CryoStor [®] cell cryopreservation medium CS10 (Sigma-Aldrich, Cat# C2874, ready-to-use) |
| Additional reagents: | <p>OptiPRO[™] SFM medium (Evercyte, Cat# MHT-033-3) supplemented with GlutaMAX[™]-I and G418 (for preparation please refer to protocol <i>Preparation of hTEC/SVTERT24-B medium</i>)</p> <p>50 µg/ml Collagen I solution (Sigma-Aldrich, Cat# C2249), for preparation please refer to protocol <i>Passaging of hTEC/SVTERT24-B cells</i></p> <p>0.05 % Trypsin-EDTA (Gibco, Cat# 25300-054, ready-to-use, stored at 4°C after thawing)</p> <p>Defined Trypsin Inhibitor (Gibco, Cat# R007100, ready-to-use, stored at 4°C after thawing)</p> <p>PBS (Sigma-Aldrich, Cat# D8537, ready-to-use, stored at RT)</p> |
| Freezing cells: | <ul style="list-style-type: none">- detach the cells from the culture vessel by using Trypsin-EDTA solution as described in protocol <i>Passaging of hTEC/SVTERT24-B cells</i>- resuspend the detached cells in growth medium and centrifuge at 170 g for 5 min- discard the supernatant- resuspend the resulting cell pellet in the remaining droplet- add freezing medium (tempered to 4°C) to reach a cell density of about 1 x 10⁶ cells/ml (for thawing in a 25 cm² culture flask)- add 1 ml of this cell suspension to each pre-cooled cryovial and immediately transfer the cells to -80°C- after 24 hours transfer the vials to the liquid nitrogen tank |
| Thawing cells: | <p>When you start cultivating the cells, please transfer the content of the original Evercyte vial containing hTEC/SVTERT24-B cells into a T25 roux flask as described in the following:</p> <ul style="list-style-type: none">- pre-coat a 25 cm² culture flask with Collagen I as described in protocol <i>Passaging of hTEC/SVTERT24-B cells</i>- add 6 ml of growth medium to the Collagen I coated 25 cm² culture flask and place the culture flask in the incubator for at least 30 min to allow the medium to reach 37°C and its normal pH- take a vial of frozen cells, rinse it outside with ethanol and pre-warm in the hand until one last piece of frozen cells is seen- then, immediately transfer the content of the vial to a 15 ml centrifugation tube filled with 9 ml of medium pre-cooled to 4°C and centrifuge for 5 min at 170 g- discard the supernatant and resuspend the cell pellet in the remaining droplet |

- add 1 ml of the pre-warmed medium to the cells, transfer the cells to the prepared culture flask and incubate at 37°C in a suitable incubator
- perform a medium change 24 hours after thawing, if the cells are already 90 % confluent at this point, they have to be passaged as described in protocol *Passaging of hTEC/SVTERT24-B cells*