

Protocol for passaging of hTEC/SVERT24-B

Version: May 2021

Evercyte Ord. No.:	CLHT-010-0024-B
Designation:	hTEC/SVERT24-B, human thymic epithelial cells
Growth medium:	<p>OptiPRO™ SFM (Gibco) supplemented with GlutaMAX™-I and G418:</p> <p><u>Final components:</u></p> <p>OptiPRO™ SFM (Gibco, Cat# 12309-019)</p> <p>2 mM GlutaMAX™-I (Gibco, Cat# 35050-038, ready-to-use, stored at RT)</p> <p>50 µg/ml G418 (InvivoGen, Cat# ant-gn5, 100 mg/ml stock solution, ready-to-use, stored at -20°C)</p> <ul style="list-style-type: none"> - take one bottle of OptiPRO™ SFM medium (500 ml) and discard 5 ml - add 5 ml of GlutaMAX™-I to the remaining 495 ml OptiPRO™ SFM medium - add 250 µl of G418 stock solution - mix properly - store at 4°C for a maximum of 4 weeks - temper the medium to RT (not 37°C) before use
Coating solution:	<p>50 µg/ml Collagen I solution</p> <p>The coating solution is prepared by mixing the following components:</p> <p>Collagen I (Sigma-Aldrich, Cat# C2249, 3 mg/ml, stored at 4°C)</p> <p>PBS (Sigma-Aldrich, Cat# D8537, stored at RT)</p> <ul style="list-style-type: none"> - take 29.5 ml PBS and transfer to a 50 ml centrifugation tube - add 0.5 ml of Collagen I solution (3 mg/ml stock solution) - mix carefully <p>For coating of a T25 cell culture flask proceed as follows:</p> <ul style="list-style-type: none"> - pipette 2 ml of diluted Collagen I solution (50 µg/ml) to a T25 roux flask - completely wet the surface of the culture flask (80 µl/cm²) - incubate the culture flask for a minimum of 30 min at 37°C - remove excess of Collagen I solution - rinse culture flask once with PBS (160 µl/cm²) - use culture flask immediately for seeding of cells, the surface must not dry out
Additional reagents:	<p>0.05 % Trypsin-EDTA (Gibco, Cat# 25300-054, stored at 4°C after thawing)</p> <p>Defined Trypsin Inhibitor (Gibco, Cat# R007100, stored at 4°C after thawing)</p>

PBS (Sigma-Aldrich, Cat# D8537, stored at RT)

Collagen I solution (Sigma-Aldrich, Cat# C2249, 3 mg/ml stock solution, stored at 4°C, prepared see above)

-
- Passaging of cells:
- remove and discard the culture medium
 - wash the cells once with PBS (each 160 $\mu\text{l}/\text{cm}^2$), remove PBS completely
 - add Trypsin-EDTA solution (20 $\mu\text{l}/\text{cm}^2$), make sure that all cells have been in contact with this solution
 - incubate the culture flask at 37°C for approximately 1-3 min
 - observe cell detachment under an inverted microscope
 - as soon as all cells are detached (if necessary, agitate the cells by gently hitting the flask), add Defined Trypsin Inhibitor (20 $\mu\text{l}/\text{cm}^2$)
 - resuspend the cells in growth medium (about 160 $\mu\text{l}/\text{cm}^2$) and aspirate the cells by pipetting
 - centrifuge at 170 g for 5 min
 - discard the supernatant, resuspend the cell pellet in the remaining droplet and add growth medium
 - transfer appropriate aliquots of the cell suspension to Collagen I pre-coated culture vessels supplemented with growth medium (final volume of 240 $\mu\text{l}/\text{cm}^2$)
 - a split ratio of 1:6 to 1:8 twice a week is recommended (after having reached about 90 %)
 - cultivate cells at 37°C in a humidified atmosphere with 5 % CO₂