

## Protocol for cryopreservation of HUVEC/TERT2

Version: September 2021

Evercyte Ord. No.:	CHT-006-0008
Designation:	HUVEC/TERT2, human umbilical vein endothelial cells
Freezing medium:	<p>EGM™ Endothelial Cell Growth Medium BulletKit™ (Lonza, Cat# CC-3124) supplemented with FBS and G418 (see protocol <i>Passaging of HUVEC/TERT2 cells</i>) 10 % DMSO (Sigma-Aldrich, Cat# D2650, ready-to-use, stored at RT)</p> <p>Preparation of 10 ml freezing medium, prepare just before use:</p> <ul style="list-style-type: none"><li>- take 9 ml of HUVEC/TERT2 growth medium and transfer to 15 ml centrifugation tube</li><li>- add 1 ml of DMSO</li><li>- mix properly and store at 4°C</li></ul>
Additional reagents:	<p>0.1 % Gelatin (Sigma-Aldrich, Cat# G1393, 2%, stored at 4°C, diluted in PBS) 0.05 % Trypsin-EDTA (Gibco, Cat# 25300-054, ready-to-use, stored at 4°C) PBS (Sigma-Aldrich, Cat# D8537, ready-to-use, stored at RT)</p>
Freezing cells:	<ul style="list-style-type: none"><li>- detach the cells from the culture vessel by using Trypsin-EDTA solution as described in protocol <i>Passaging of HUVEC/TERT2 cells</i></li><li>- resuspend the detached cells in growth medium and centrifuge at 170 g for 5 min</li><li>- discard the supernatant</li><li>- resuspend the resulting cell pellet in the remaining droplet</li><li>- add freezing medium (tempered to 4°C) to reach a cell density of about <math>5 \times 10^5</math> cells/ml (for thawing in a 25 cm<sup>2</sup> culture flask)</li><li>- add 1 ml of this cell suspension to each pre-cooled cryovial and immediately transfer the cells to -80°C</li><li>- after 24 hours transfer the vials to the liquid nitrogen tank</li></ul>
Thawing cells:	<p><b>When you start cultivating the cells, please transfer the content of the original Evercyte vial containing HUVEC/TERT2 cells into a T25 roux flask as described in the following:</b></p> <ul style="list-style-type: none"><li>- pre-coat a 25 cm<sup>2</sup> culture flask with Gelatin solution as described in protocol <i>Passaging of HUVEC/TERT2 cells</i></li><li>- add 6 ml of growth medium to a 25 cm<sup>2</sup> culture flask and place the culture flask in the incubator for at least 30 min to allow the medium to reach 37°C and its normal pH</li><li>- take a vial of frozen cells, rinse it outside with Ethanol and pre-warm in the hand until one last piece of frozen cells is seen</li><li>- then, immediately transfer the content of the vial to a 15 ml centrifugation tube pre-filled with 9 ml of medium pre-cooled to 4°C and centrifuge for 5 min at 170 g</li><li>- discard the supernatant and resuspend the cell pellet in the remaining droplet</li><li>- add 1 ml of the pre-warmed medium to the cells, transfer the cells to the prepared culture flask and incubate at 37°C in a suitable incubator</li></ul>

- perform a medium change 24 hours after thawing, if the cells are already near confluent at this point, they have to be passaged as described in protocol *Passaging of HUVEC/TERT2 cells*

---

Related products: HUVEC/TERT66, umbilical vein endothelial cells (Evercyte, Cat# CHT-006-0066)  
HDMVEC/TERT164-B, dermal microvascular endothelial cells, lymphatic (Evercyte, Cat# CHT-013-0164-B)