

Protocol for passaging of HUVEC/TERT2

Version: September 2021

| | |
|--------------------|--|
| Evercyte Ord. No.: | CHT-006-0008 |
| Designation: | HUVEC/TERT2, human umbilical vein endothelial cells |
| Growth medium: | EGM TM Endothelial Cell Growth Medium BulletKit TM (Lonza, Cat# CC-3124) supplemented with FBS and G418: |

Final components:

EBMTM basal medium (Lonza, Cat# CC-3121)
 Components of EGMTM SingleQuotTM Supplements (Lonza, Cat# CC-4133: BBE, hEGF, Hydrocortisone, Ascorbic acid)
 10 % FBS (PAN-Biotech, Cat# P30-3031)
 20 µg/ml G418 (InvivoGen, Cat# ant-gn-5, 100 mg/ml stock solution, ready-to-use)

- take one bottle of EBM basal medium (500 ml)
- discard 50 ml of EBM basal medium
- add 500 µl of Hydrocortisone (EGMTM SingleQuotTM Supplements)
- add 500 µl of hEGF (EGMTM SingleQuotTM Supplements)
- add 2 ml of BBE (EGMTM SingleQuotTM Supplements)
- add 500 µl of Ascorbic Acid (EGMTM SingleQuotTM Supplements)
- add 50 ml of FBS
- add 100 µl of G418 stock solution

- mix properly
- store at 4°C for a maximum of 4 weeks
- temper the medium to room temperature (not 37°C) before use

| | |
|----------|---|
| Coating: | 0.1 % Gelatin solution |
| | <p>The coating solution is prepared by mixing the following components: Gelatin (Sigma-Aldrich, Cat# G1393, 2 %, stored at 4°C) PBS (Sigma-Aldrich, Cat# D8537, ready-to-use, stored at RT)</p> <ul style="list-style-type: none"> - liquify the Gelatin stock solution (2 % - stored in 2 ml aliquots at 4°C) at 37°C - add 38 ml of PBS to 2 ml Gelatin stock solution (2 %) - mix carefully, store diluted Gelatin solution (0.1 %) at 37°C <p>For coating of a T25 roux flask proceed as follows:</p> <ul style="list-style-type: none"> - transfer 2 ml of Gelatin solution (0.1 %) to a T25 roux flask (final 80 µl/cm²) |

- completely wet the surface of the culture flask
- incubate the culture flask at 37°C between 10-60 min
- remove excess of Gelatin solution
- use culture flask immediately for seeding of cells, the surface must not dry out

Additional reagents: 0.1 % Gelatin (Sigma-Aldrich, Cat# G1393, 2 %), diluted in PBS
PBS (Sigma-Aldrich, Cat# D8537)
0.05 % Trypsin-EDTA (Gibco, Cat# 25300-054)

- Passaging of cells:
- the new culture flasks have to be coated as described above
 - remove and discard the culture medium
 - wash the cells once with PBS (each 160 µl/cm²), remove PBS completely
 - add Trypsin-EDTA solution (20 µl/cm²), make sure that all cells have been in contact with this solution
 - incubate the culture flask at 37°C for approximately 3 min
 - observe cell detachment under an inverted microscope
 - as soon as all cells are detached (if necessary, agitate the cells by gently hitting the flask), add growth medium (about 160 µl/cm²)
 - add appropriate aliquots of the cell suspension to Gelatin pre-coated culture vessels supplemented with growth medium (final volume of 240 µl/cm²)
 - a split ratio of 1:8 twice a week is recommended (after having reached about 90-95 % confluence)
 - cultivate cells at 37°C in a humidified atmosphere with 5% CO₂

Related products: HUVEC/TERT66, umbilical vein endothelial cells (Evercyte, Cat# CHT-006-0066)
HDMVEC/TERT164-B, dermal microvascular endothelial cells, lymphatic (Evercyte, Cat# CHT-013-0164-B)
