

Protocol for passing of NHEK/SVTERT3-5

page 1 of 2

Version: February 2024

Evercyte Ord. No.:	CLHT-011-0026-5
Designation:	NHEK/SVTERT3-5, human epidermal keratinocytes
Growth medium:	<p>KGM™ Gold Keratinocyte Growth Medium BulletKit™ (Lonza, Cat# 00192060) supplemented with G418:</p> <p><u>Final components:</u> KBM™ Gold™ Basal Medium (Lonza, Cat# 00192151) Components of KGM™ Gold™ SingleQuots™ supplements (Lonza, Cat# 00192152: BPE, hEGF, Insulin, Hydrocortisone, Transferrin, Epinephrine, without GA) 50 µg/ml G418 (InvivoGen, Cat# ant-gn5, 100 mg/ml stock solution, ready-to-use)</p> <ul style="list-style-type: none">- take one bottle of KBM™ Gold™ Basal Medium (500 ml)- add 2 ml of BPE (component of KGM™ Gold™ SingleQuots™)- add 500 µl of hEGF (component of KGM™ Gold™ SingleQuots™)- add 500 µl of Insulin (component of KGM™ Gold™ SingleQuots™)- add 500 µl of Hydrocortisone (component of KGM™ Gold™ SingleQuots™)- add 500 µl of Transferrin (component of KGM™ Gold™ SingleQuots™)- add 250 µl of Epinephrine (component of KGM™ Gold™ SingleQuots™)- add 250 µl of G418 stock solution (100 mg/ml)- mix properly- store at 4°C for a maximum of 4 weeks- temper the medium to room temperature (not 37°C) before use
Additional reagents:	<p>0.05% Trypsin-EDTA (Gibco, Cat# 25300-054, ready-to-use, stored at 4°C) Defined Trypsin Inhibitor (Gibco, Cat# R007100, ready-to-use, stored at 4°C) PBS (Sigma-Aldrich, Cat# D8537, ready-to-use, stored at RT)</p>
Passaging of cells:	<ul style="list-style-type: none">- remove and discard the culture medium- wash the cells once with PBS (160 µl/cm²), remove PBS completely- add Trypsin-EDTA solution (20 µl/cm²), make sure that all cells have been in contact with this solution- incubate the culture flask at 37°C for approximately 2-3 min, cells have typically rounded up but have not detached yet- tap side of the flask to apply mechanical force and further incubate for 2-3 min- observe cell detachment under an inverted microscope- as soon as all cells are detached (if necessary, agitate the cells by gently hitting the flask), add Defined Trypsin Inhibitor (20 µl/cm²)

- resuspend the cells in growth medium (about 160 $\mu\text{l}/\text{cm}^2$) and aspirate the cells by pipetting
 - pass cells through EASY strainer 100 μm (Greiner Cat# 542000)
 - centrifuge at 170 g for 5 min
 - discard the supernatant, resuspend the cell pellet in the remaining droplet and add growth medium
 - transfer appropriate aliquots of the cell suspension to new culture vessels supplemented with growth medium (final volume of 240 $\mu\text{l}/\text{cm}^2$)
 - a split ratio of 1:4 to 1:6 twice a week is recommended (after having reached about 60-70 % confluence)
 - perform a medium change after 2-3 days if cells have not reached required cell density
 - cultivate cells at 37°C in a humidified atmosphere with 5% CO_2
-