

## Protocol for preparation of EN-MSC/TERT352-B medium

page 1 of 1

Version: March 2025

Evercyte Ord. No.: CHT-060-0352-B

---

Designation: EN-MSC/TERT352-B, human endometrium derived mesenchymal stromal cells

---

Growth medium: RoosterNourish™-MSC-XF (RoosterBio, Cat# K82016) supplemented with G418:

Final components:

RoosterBasal™2.0-CC (RoosterBio, Cat# M22520)

RoosterBooster™-MSC-XF (RoosterBio, Cat# SU-016)

200 µg/ml G418 (InvivoGen, Cat# ant-gn5, ready-to-use, 100 mg/ml)

- take one bottle of RoosterBasal™2.0-CC medium (500 ml)
  - add 10 mL RoosterBooster™-MSC-XF
  - add 1 ml of G418 stock solution
  
  - mix properly
  - store at 4°C for a maximum of 2 weeks
  - temper the medium to room temperature (not 37°C) before use
-